

proVi⁺



PRODUCTS
SYRINGES

«COMMITTED TO QUALITY SINCE 2014»

med inject

medInject Syringes are manufactured to achieve the highest level of accuracy and precision possible. The needles are precisely honed and beveled for ease of penetration and less patient trauma. The bold scale markings on the barrel assure the precise measurement of the medication dose. An integrated plunger stop helps prevent the plunger from being accidentally pulled out. The ergonomically designed finger grips and push-button allow a perfectly controlled and stable injection.



LUER LOCK SYRINGES

- Sterile and Latex free
- Clear scale markings to enable precise dosage
- Colored plunger for better visualization

Available sizes

3ml - 23G (1") / 23G (1 1/4"), 5ml - 23G (1") / 23G (1 1/4")
10ml 21G (1 1/2"), 20ml 21G (1 1/2"), 50ml flushing syringe
without needle - purpose of irrigation

med inject

medInject Syringes for insulin injection provide the highest quality from plunger to tip. With 31G needles we want to minimize trauma and discomfort, which is especially important when the patient needs injections on a daily basis. The permanently attached needle limits dead space. This results in less waste of costly insulin and facilitates air bubble removal. Every needle is individually inspected to ensure the patient experience a more comfortable injection.



INSULIN SYRINGES

- Permanently attached needle
- Beveled tip to minimize discomfort
- Smooth plunger action for greater filling accuracy

Available sizes

0.5ml 31G x 5mm, 1ml 31G x 5mm

med comfort

The sterile and single-use **medComfort Insulin Pen Needles** are intended to be used with a pen injector device for the subcutaneous injection of insulin. They are commonly used by people with diabetes who often require multiple daily insulin injections. The pen needles come in a variety of lengths and widths to suit all body types and are designed for a comfortable injection experience supporting patient satisfaction.



INSULIN

PEN NEEDLES

- Single use
- Non toxic
- Non pyrogenic

Available sizes

32G x 4mm, 31G x 6mm, 31G x 8mm

med connect

The **medConnect Intravenous Cannulas** are inserted into a vein, primarily for the administration of intravenous fluids, for obtaining blood samples and for administering medicines. Two wings are ensuring easy and secure positioning on the skin of the patient. medConnect Venous Cannula are available in all common gauge sizes. Injection port is equipped with a one way valve with a flange attached to the catheter/shaft. Injection port caps are color coded for different sizes.



INTROVENOUS

CANNULAS

- Lorem Ipsum
- Lorem Ipsum
- Lorem Ipsum

Available sizes

32G x 4mm, 31G x 6mm, 31G x 8mm

proVi⁺



ProVi Swiss GmbH

Widnauerstrasse 1
9435 Heerbrugg
Switzerland

+41 71 544 43 24
info@proviswiss.com
www.proviswiss.com