

MONOFILAMENT POLYAMIDE SUTURE

NONABSORBABLE SURGICAL SUTURES













Nylon sutures are non-absorbable sutures and possess excellent tensile strength. Nylon sutures are available in black colour. Nylon sutures have excellent knot security properties and can be easily removed with no tissue adherence. These sutures are infection resistant.











Nylon sutures are widely used for general closure, skin and plastic surgery. They are usually not recommended for attachment of artificial prostheses in cardiovascular surgery. They also do not support infection and maintain tensile strength indefinitely in tissues.

Virtues	Being monofilament, Linex is remarkably smooth and gives excellent knot security. Possesses good tensile strength. Easy removal with no tissue adherence.
Material	Polyamide 6 suture (Nylon). Blue/Black coloured for easy identification in the operating field.
Tissue Reaction	Grade-0 Reactivity - None
Applications	General closure, Skin and Plastic Surgery. Also used for Ophthalmic and Micro Surgeries.
Choice	Loop, Needled and pre-cut sterile sutures
Caution	Not recommended for attachment of artificial prostheses in Cardiovascular surgery.













NONABSORBABLE SURGICAL SUTURES

NEEDLE LENGTH	NEEDLE DESCRIPTION	CUT	SUTURE LENGTH	SUTURE SIZE					
				10-0(0.2)	6-0(0.7)	5-0(1)	4-0(1.5)	3-0(2)	2-0(3)
 6 mm	3/8 Circle spatulated micro point		38 cm	SC 3719					
 10 mm	3/8 Circle Reverse Cutting		70 cm		SC 3320				
 12 mm	3/8 Circle Reverse Cutting		70 cm			SC 3317			
 16 mm	3/8 Circle Cutting		35 cm			SC 3315			
			70 cm				SC 3318	SC 3321	
 26 mm	3/8 Circle Reverse Cutting		35 cm					SC 3328S	SC 3329S
			70 cm					SC 3328M	SC 3329M
			100 cm					SC 3328	SC 3329
 30 mm	1/2 Circle Round Body		100 cm						SC 3720

NEEDLE LENGTH	NEEDLE DESCRIPTION	CUT	SUTURE LENGTH	SUTURE SIZE		
				2-0(3)	0(3.5)	1(4)
 36 mm	3/8 Circle Cutting		100 cm	SC 3236	SC 3239	
 40 mm	1/2 Circle Reverse Cutting		100 cm	SC 3335	SC 3337	
			45 cm	SC 3335S	SC 3337S	
 40 mm	1/2 Circle Round Body		150 cm Loop		SC 3340	
 40 mm Heavy	1/2 Circle Reverse Cutting		100 cm		SC 3339	
 40 mm Heavy	1/2 Circle Round Body		100 cm		SC 3346	SC 3347

NONABSORBABLE SURGICAL SUTURES

NEEDLE LENGTH	NEEDLE DESCRIPTION	CUT	SUTURE LENGTH	SUTURE SIZE		
				2-0(3)	0(3.5)	1(4)
 45 mm	3/8 Circle Reverse Cutting		45 cm	SC 3336S		
			70 cm	SC 3336M		
			100 cm	SC 3336		SC 3334
 45 mm Heavy	1/2 Circle Reverse Cutting		45cm			SC 3338S
			100 cm			SC 3338
 50 mm Heavy	1/2 Circle Round Body		150 cm Loop			SC 3348
 60 mm	Straight Cutting		70 cm	SC 3390	SC 3389	
 90 mm	3/8 Circle Cutting		100 cm			SC 3397

PRE-CUT SUTURES

SUTURE LENGTH	PACKING	SUTURE SIZE		
		2-0(3)	0(3.5)	1(4)
5 x 35 cm	Sterilised Monofilament Polyamide in sterilised foil over wrap. 5 x 35 cm box of 12's pack	SC 703	SC 704	SC 705
2 x 76 cm	Sterilised Monofilament Polyamide in sterilised foil over wrap. 5 x 35 cm box of 12's pack	SC 903	SC 904	SC 905