

BLACK BRAIDED SILK SUTURE

NONABSORBABLE SURGICAL SUTURES

















Silk sutures are entwined thread from silk worm larva. Silk sutures are coated and braided non-absorbable sutures. Silk sutures are specially treated with silicone to give strength and smoothen the surface which provides easy passage through tissues.















Silk sutures are widely used as ligature and are also used for other applications like skin, ophthalmic, GI tract, etc. and are not indicated for use in biliary and urinary tract surgery.

Virtues	Most widely used. Excellent handling characteristics. Good knot security.
Material	Entwined thread from silk worm larva. Treated with silicone to give strength and smooth surface.
Tissue Reaction	Grade-0 Reactivity - None
Applications	General closure, G. I. Tract, Plastic, Skin, Ophthalmic, Cardiovascular Surgeries. Widely used as ligature.
Choice	Needled and pre-cut sterile sutures
Caution	Avoid using in the presence of infection. Not used in biliary and urinary tract surgery.











NONABSORBABLE SURGICAL SUTURES

NEEDLE LENGTH	NEEDLE DESCRIPTION	CUT	SUTURE LENGTH	SUTURE SIZE			
				5-0(1)	4-0(1.5)	3-0(2)	2-0(3)
 12 mm	3/8 Circle Reverse Cutting		76 cm	SCR 5080			
 16 mm	3/8 Circle Cutting		76 cm		SCR 5001	SCR 5002	
 16 mm	3/8 Circle Round Body		76 cm	SCR 5081	SCR 5082		
 16 mm	1/2 Circle Round Body		76 cm		SCR 5049		
 20 mm	1/2 Circle Round Body		76 cm		SCR 5086	SCR 5087	
 25 mm	1/2 Circle Round Body		76 cm			SCR 5070	
 25 mm	1/2 Circle Taper Cutting		76 cm				SCR 5670
 26 mm	3/8 Circle Reverse Cutting		76 cm			SCR 5028	SCR 5029

NEEDLE LENGTH	NEEDLE DESCRIPTION	CUT	SUTURE LENGTH	SUTURE SIZE			
				2-0(3)	0(3.5)	1(4)	2(5)
 30 mm Heavy	1/2 Circle Cutting		76 cm	SCR 5072			
 30 mm	1/2 Circle Cutting		76 cm		SCR 5077		
 30 mm	1/2 Circle Round Body		76 cm	SCR 5333	SCR 5334		
			90 cm	SCR 5331			
 36 mm	3/8 Circle Cutting		76 cm		SCR 5360	SCR 5361	SCR 5362
 40 mm	1/2 Circle Cutting		76 cm		SCR 5064		
 40 mm	1/2 Circle Reverse Cutting		76 cm	SCR 5067	SCR 5068	SCR 5069	
 40 mm	1/2 Circle Round Body		76 cm		SCR 5075		SCR 5073

NONABSORBABLE SURGICAL SUTURES

NEEDLE LENGTH	NEEDLE DESCRIPTION	CUT	SUTURE LENGTH	SUTURE SIZE		
				2-0(3)	0(3.5)	1(4)
 40 mm Heavy	1/2 Circle Round Body		90 cm			SCR 5063
 45 mm	3/8 Circle Reverse Cutting		90 cm	SCR 5036	SCR 5037	
 50 mm	1/2 Circle Reverse Cutting		76 cm	SCR 5052		
 60 mm	3/8 Circle Cutting		76 cm			SCR 5062

PRECUT SUTURES (NON NEEDED)

SUTURE LENGTH	PACKING	SUTURE SIZE				
		3-0(2)	2-0(3)	0(3.5)	1(4)	2(5)
2 x 76 cm	Sterilised black braided silk in sterile foil over wrap. 2 x 76 cm box of 12's pack.	SCR 212	SCR 213	SCR 214	SCR 215	SCR 216

SILK REELS

SUTURE LENGTH 25 METERS	PACKING	SUTURE SIZE				
		3-0(2)	2-0(3)	0(3.5)	1(4)	2(5)
6 Reels in a Box Non sterile	Surgical black braided 6 reels x 25 meters. (Non-sterile)	SCR 822	SCR 823	SCR 824	SCR 825	SCR 826